

## Training Blood Collection Pad



This life-like blood collection pad provides opportunity to practice blood sampling and venous catheterisation in small animals.

**Product codes:**

Reference: VTQ-INJ-PAD-SA

**Product description:**

This blood collection pad has two integrated veins, perfectly mimicking the characteristics of subcutaneous blood vessels.

The two veins simulate different levels of blood stasis and thus allow training at various levels of difficulty, as is often the case in real-life conditions.

Our suture pads are reusable, tear-resistant and retain their shape and texture even after repeated use.

Our practical pad holder can be used to mount the suture pads in a vertical or a horizontal position, as required.

**Product features:**

Dimensions: 15 x 16 cm

## Available models

