

## Rodent thermometer BIO-TK8851



The thermometer TK8851 is combining accuracy of the digital device and simplicity of reading with a large readable screen - Entry level tool for all your temperature measurements on rodents (rats or mice).

Probes for rat or mice, of your choice : not included in the thermometer price

### Product codes:

Reference: BIO-TK8851

### Product description:

**Thermometer BIO-TK8851** is specially targeting measurements on small rodents like rats and mice, and offers a wide area of functionalities making it the ideal thermometer for laboratory applications during your studies involving animal models.

The T1 channel offers a resolution of 0.1°C, while the large integrated LCD display is able to display 2.5 measurements/second. With its advanced functionalities (records max, records min, hold), the TK8851 thermometer for rodents is an **easy-to-use instrument yet offering state-of-the art performance for your research.**

Its solid construction and rubber keys, along with the plastic case, make the TK8851 a **very durable instrument** fit for all applications and conditions.

### Product features:

Weight: 150g

Battery: 9V

Accuracy: Temperature : +/- (0.1% reading + 0.7°C) (for 18 to 28°C ambient temp.)

Temperature distribution: Operating : 0 to 50°C

Temperature range: -200°C to 390°C (according to probe) (for 18 to 28°C ambient temp.)

### Available models

