

Spectrafuge 24D Digital Microcentrifuge

The high end of the performance spectrum



Product codes:

Reference: C-2400

Product description:

An easy-access 24 x 1.5 ml rotor is supplied with the Spectrafuge 24D.

This innovative rotor places the tubes on a plateau, where the user can easily grip their tops. The powerful brushless motor quickly accelerates the rotor to a maximum speed of 13,300rpm/16,000 x g.

Product features:

Dimensions: 23.5 x 29.3 x 21.6 cm

Time: 0 to 30 minutes or continuous "Quick" button for momentary operation

Speed: 500rpm to 13,300rpm

Capacity: 24 x 1.5/2.0 ml

Centrifugal force: 16,300 x g

Available models

Color: Blue (C-2400-B), Teal (C-2400-T), Purple (C-2400-P), Red (C-2400-R)

O LAB