

## Renasil Silicone Rubber Tubing



An essentially non reactive platinum cured silicon tubing

**Product codes:**

Reference: SIL

**Product description:**

This high grade silicone rubber tubing is now available to the research community.

**Product features:**

Length: 7.60 m

**Available models**

Diameter: 0.64 OD x 0.30 ID mm (SIL025), 0.94 OD x 0.64 ID mm (SIL037), 1.19 OD x 0.64 ID mm (SIL047), 1.65 OD x 0.76 ID mm (SIL065), 2.03 OD x 1.02 ID mm (SIL-080)