


# PTFE Heat Shrink Tubing

USP Class VI Approved

**Product codes:**

Reference: PTFE



## Product description:

Medical grade Polytetrafluoroethylene tubing has high tensile strength and preserves high fidelity pressure waveforms with very small internal diameter tubing.

## Product features:

Length: 30.5 cm

## Available models

Size: Exp. ID 0.97 mm - Rec. ID 0.46 mm - Rec. wall 0.23 mm - Tol. +/- 0.051 mm (PTFE-018),  
Exp. ID 1.5 mm - Rec. ID 0.99 mm - Rec. wall 0.3 mm - Tol. +/- 0.076 mm (PTFE-039)