


Universal Restrainer



Experiments Made Easier - different models for rats and mice, with smoked or clear canvas

Product codes:

Reference: 500M-C

Product description:

The Universal Restrainer was designed to allow the handler to perform all sampling and injection procedures in a single immobilization.

- the slit in the abdominal ventral area of the rat allows : access to the bladder, gastrointestinal catheters, fistulas, intraperitoneal administration, cardiac punctures
- the dorsal slit allows : the insertion of catheters, subcutaneous injections, EKG and EEG
- the slit at the tail allows : for sampling and injection, heart puncture

Made of Stainless steel and plastic

Product features:

Weight: rat up to 350g

Weight: mice 20g to 60g

Available models

Color: Clear (500M-C), Smoked (500M-S)

Animal: Mice 20g to 60g (500M-C), Mice 20g to 60g (500M-S)

