

Ideal Micro-Drill

Designed for research applications that require surgical burrs and trephines

**Product codes:**

Reference: MD-1200

Product description:

Designed for research applications that require surgical burrs and trephines. The drill is constructed of light-weight aluminum alloy for balance and control and is powered by a rechargeable 6-volt nickel cadmium battery. Under normal operating conditions the unit will function for 8 hours between charges.

Product features:

Power supply: 6-Volts DC
Length: 17.5 cm (Without Burr)
Diameter: 1.9 mm
Speed: 12,000 rpm

Available models