

Re-usable Feeding Needles



A wide range of re-usable feeding needles

Product codes:

Reference: N-

Product description:

Re-usable feeding needle.

Stainless steel

Select the proper feeding needle for your specific need

ANIMAL/ WEIGHT RANGE/ BALL DIAMETER

Mice 25 - 125 g 1.25 mm

Rats 100 - 400 g 2.25 mm

Rats 350 - 500 g 3.00 mm

Other over 450 g 4.00 mm

Product features:

Available models

Size: 24ga - 25.4mm Straight - Ball dia. 1.25mm (N-PK-001), 22ga - 25.4mm Straight - Ball dia. 1.25mm (N-PK-002), 20ga - 76.2mm Straight - Ball dia. 2.25mm (N-PK-003), 18ga - 76.2mm Straight - Ball dia. 2.25mm (N-PK-004), 20ga - 76.2mm Curved - Ball dia. 2.25mm (N-PK-005), 18ga - 76.2mm Curved - Ball dia. 2.25mm (N-PK-006), 16ga - 76.2mm Curved - Ball dia. 3mm (N-PK-007), 18ga - 25.4mm Straight - Ball dia. 2.25mm (N-PK-008), 16ga - 101.6mm Curved - Ball dia. 3mm (N-PK-009), 20ga - 38.1mm Curved - Ball dia. 2.25mm (N-PK-010), 22ga - 38.1mm Straight - Ball dia. 1.25mm (N-PK-020)

Units: 4 per pack (N-PK-001), 4 per pack (N-PK-002), 4 per pack (N-PK-003), 4 per pack (N-PK-004), 4 per pack (N-PK-005), 4 per pack (N-PK-006), 4 per pack (N-PK-007), 4 per pack (N-PK-008), 4 per pack (N-PK-009), 4 per pack (N-PK-010), 4 per pack (N-PK-020)

