

## Tube for Tail Illuminator



Tube only

**Product codes:**

Reference: MTI-TUBE1

### Product description:

Spare Tube for MTI  
Easy to clean

**Ref Int. Diameter  
in cm Int. Diameter**

**in "** MTI-TUBE1 3,18 1.25 " - Standard MTI-TUBE2 2,86 1.12 " - Small MTI-TUBE3 rotating MTI-TUBE4 3,81 1.5 " - Large

### Product features:

### Available models

Diameter: 3.18 cm (MTI-TUBE1), 2.86 cm (MTI-TUBE2), 3.81 cm (MTI-TUBE4), rotating (MTI-TUBE3)