


Rat Jugular Vein Catheter



New Polyurethane (Micro-Renathane and RenaPulse) Rat Vascular Catheters. 3 series are available: Rat Jugular Vein Catheter, Rat Thoracic Carotid Artery and Rat Femoral Artery Catheter.

Product codes:

Reference: R-JVC-

Product description:

Total catheter length: 120 mm
Intravascular length: 32 mm
Port length: 25 mm
Number of retention beads: 3

Product features:**Available models**

Composition: MRE-037 (R-JVC-M37), MRE-037 - Tip rounding (R-JVC-M37-R)