


Polyethylene tubing



Ideal for short-term catheterizations and for ex vivo applications

Product codes:

Reference: PE

Product description:**Product features:****Available models**

Length: 50 feet (15.24 m) (PE-10-CL-50), 50 feet (15.24 m) (PE-20-CL-50), 50 feet (15.24 m) (PE-30-CL-50), 50 feet (15.24 m) (PE-50-CL-50), 50 feet (15.24 m) (PE-60-CL-50), 50 feet (15.24 m) (PE-100-CL-50), 100 feet (30.5 m) (PE-10-CL-100), 100 feet (30.5 m) (PE-20-CL-100), 100 feet (30.5 m) (PE-30-CL-100), 100 feet (30.5 m) (PE-50-CL-100), 100 feet (30.5 m) (PE-60-CL-100), 100 feet (30.5 m) (PE-100-CL-100), 200 feet (61 m) (PE-10-CL-200), 200 feet (61 m) (PE-20-CL-200), 200 feet (61 m) (PE-30-CL-200), 200 feet (61 m) (PE-50-CL-200), 200 feet (61 m) (PE-60-CL-200), 200 feet (61 m) (PE-100-CL-200), 500 feet (152.5 m) (PE-10-CL-500), 500 feet (152.5 m) (PE-20-CL-500), 500 feet (152.5 m) (PE-30-CL-500), 500 feet (152.5 m) (PE-50-CL-500), 500 feet (152.5 m) (PE-60-CL-500), 500 feet (152.5 m) (PE-100-CL-500), 1000 feet (305 m) (PE-10-CL-1000), 1000 feet (305 m) (PE-20-CL-1000), 1000 feet (305 m) (PE-30-CL-1000), 1000 feet (305 m) (PE-50-CL-1000), 1000 feet (305 m) (PE-60-CL-1000), 1000 feet (305 m) (PE-100-CL-1000)

Diameter: 0.60 OD x 0.28 mm ID (PE-10-CL-50), 1.09 OD x 0.38 mm ID (PE-20-CL-50), 0.76 OD x 0.43 mm ID (PE-30-CL-50), 0.965 OD x 0.58 mm ID (PE-50-CL-50), 1.22 OD x 0.76 mm ID (PE-60-CL-50), 1.52 OD x 0.86 mm ID (PE-100-CL-50), 0.60 OD x 0.28 mm ID (PE-10-CL-100), 1.09 OD x 0.38 mm ID (PE-20-CL-100), 0.76 OD x 0.43 mm ID (PE-30-CL-100), 0.965 OD x 0.58 mm ID (PE-50-CL-100), 1.22 OD x 0.76 mm ID (PE-60-CL-100), 1.52 OD x 0.86 mm ID (PE-100-CL-100), 0.60 OD x 0.28 mm ID (PE-10-CL-200), 1.09 OD x 0.38 mm ID (PE-20-CL-200), 0.76 OD x 0.43 mm ID (PE-30-CL-200), 0.965 OD x 0.58 mm ID (PE-50-CL-200), 1.22 OD x 0.76 mm ID (PE-60-CL-200), 1.52 OD x 0.86 mm ID (PE-100-CL-200), 0.60 OD x 0.28 mm ID (PE-10-CL-500), 1.09 OD x 0.38 mm ID (PE-20-CL-500), 0.76 OD x 0.43 mm ID (PE-30-CL-500), 0.965 OD x 0.58 mm ID (PE-50-CL-500), 1.22 OD x 0.76 mm ID (PE-60-CL-500), 1.52 OD x 0.86 mm ID (PE-100-CL-500), 0.60 OD x 0.28 mm ID (PE-10-CL-1000), 1.09 OD x 0.38 mm ID (PE-20-CL-1000), 0.76 OD x 0.43 mm ID (PE-30-CL-1000), 0.965 OD x 0.58 mm ID (PE-50-CL-1000), 1.22 OD x 0.76 mm ID (PE-60-CL-1000), 1.52 OD x 0.86 mm ID (PE-100-CL-1000)

LAB