

Self-regulated heating pad

The best choice if you need a self-regulated heating pad with a tight budget

**Product codes:**

Reference: BIO-HE002

Product description:

Our heating pad can be preset within the range 10°C to 50°C.

Product features:

Weight: 3 kg
Dimensions: 55cm x 33cm and 37cm x 37cm
Accuracy: +/- 0,5°C
Heating up time: Approx. 10 min
Power supply: 220V
Temperature range: 10 °C - 50 °C
Resolution: 0.1 °C

Available models

Size: L : 55cm x 33cm (BIO-HE002), M : 37cm x 37cm (BIO-HE002M)