


Rad Disk and Rad Disk Mini



Polycarbonate rodent microisolation irradiator cages

Product codes:

Reference: IRD-P

Product description:

Because Rad disk and Rad disk mini are made entirely of polycarbonate plastic, they are autoclavable. The two cage sizes are designed to fit into radiation chambers 5" high and 9" diameter and larger. Used as instructed, Rad disk and Rad disk mini prevent contamination of mice before, during and after irradiation. They are also a versatile carrier to protect rodents while transporting them in and out of facilities.

Cage is sold with a filter

Product features:

Available models

Size: 8.875" D x 3.875" H (22.5 cm D x 9.8 cm H) (IRD-P-M), 10.75" D x 4" H (27.3 cm D x 10.2 cm H) (IRD-P-R)