


# Red Short Base (Shorti) Tail Vein Injection Tube for mice



Tail vein restraint system, mounted on a wider base plate allowing the tail to drop vertically for greater range of motion for venipuncture. For mice

**Product codes:**

Reference: SHORTI-STD

**Product description:**

Tail vein restraint system, mounted on a wider baseplate **allowing the tail to drop vertically for greater range of motion for venipuncture.**

**Product features:**

Dimensions: Base : 9.4 x 4.3 cm

Length: Tube : 8.6 cm

Diameter: Tube : 3.18 cm

**Available models**