

## Mouse restrainer for BP-2000



to be used on BP-2000 unit only

**Product codes:**

Reference: BP-MH

### **Product description:**

Restrainer for mice, for BP-2000

Feature a magnet base to keep the mouse at the correct place on the tray (tray not included)

Black color and holes at the front side keep the animal quiet for the duration of the blood pressure test.

Easy to clean and urine-resistant

### **Product features:**

Material: anodized metal

Colors: black

### **Available models**

Animal Size: very small mice (BP-MH0), mice under 30 grams (BP-MH1), mice 30-45 grams (BP-MH2), mice larger than 45 grams (BP-MH3)

