


TCAT Animal Heating Plate



Maintain temperature during surgery

Product codes:

Reference: HP-

Product description:

Current heating plate for small rodents (mice and rats), used to maintain animal body temperature during surgery procedures.

Reference Animal Dimensions (cm) HP-1M Rat 11,1 W x 0,6 H x 21,6 L HP-4M Mice 3,8 W x 0,6 H x 10,2 L

Standard plates are **black anodized aluminum**

Both models have small feet attached giving them an elevation of 0,6 cm above the surface it sits on, such as a table or stereotaxic stage.

They are easily cleaned with alcohol.

They are provided with a **feedback sensor to prevent overheating when used with our TCAT-2 series of controllers.**

Product features:

Available models

Animal: Rat (HP-1M), Mouse (HP-4M)

