

Sutures with needles 12/pk



Silk suture spool with swagged needles
pack of 12 needles

Product codes:

Reference: SUT

Product description:

SUT-1074: 3/8 circle reverse cutting, 19mm needle

3/0 silk black braided, length : 45,5cm (18")

SUT-7116: 1/2 circle taper point , 17,4mm needle

4/0, length 60,9 cm (24")

SUT-9403: Straight taper, slight angled needle, 44,4mm

3/0, length : 76,2cm (30")

Product features:

Material: Silk

Available models

Size: USP 3/0 - 45,5cm (SUT-1074), USP 4/0 - 61cm (SUT-7116)

