

## La Pias Ear Tag - One color



La Pias Ear tag, box of 100, numbered - single color of your choice

**Product codes:**

Reference: EP-56780

**Product description:**

Specially suited for use with small lab and research animals (mice, rats and other rodents). The light weight easy-to-read tags allow for trouble-free identification.

- Ear tags are **numbered from 1-100 (available up to 999)**
- **available in 5 different colors : Green, Blue, Red, Light blue, Black**
- also available as Stainless Steel (non-colored)
- Once applied, they cannot be easily removed
- Laser-etched number provides durability
- ear tag are autoclavable (applicator not)

**Packing : one box of 100 pce****Material**

Numbered Tag: Aluminum  
Clasp: Stainless Steel

**Size**

Mouse: 3.5mm O.D. with 3mm Pin

Rat: 3.5mm O.D. with 5mm Pin

**Product features:**

Material: Tag : Aluminium / Clasp : Stainless steel

**Available models**

Color: Red (EP-56780-100), Blue (EP-56782-100), Teal (EP-56783-100), Green (EP-56781-100), Black (EP-56784-100), Red (EP-56780-200), Blue (EP-56782-200), Teal (EP-56783-200), Green (EP-56781-200), Black (EP-56784-200), Red (EP-56786-100), Blue (EP-56788-100), Teal (EP-56789-100), Black (EP-56790-100), Green (EP-56787-100), Red (EP-56786-200), Blue (EP-56788-200), Teal (EP-56789-200), Black (EP-56790-200), Green (EP-56787-200)

Animal: Mouse (EP-56780-100), Mouse (EP-56782-100), Mouse (EP-56783-100), Mouse (EP-56781-100), Mouse (EP-56784-100), Mouse (EP-56780-200), Mouse (EP-56782-200), Mouse (EP-56783-200), Mouse (EP-56781-200), Mouse (EP-56784-200), Rat (EP-56786-100), Rat (EP-56788-100), Rat (EP-56789-100), Rat (EP-56790-100), Rat (EP-56787-100), Rat (EP-56786-200), Rat (EP-56788-200), Rat (EP-56789-200), Rat (EP-56790-200), Rat (EP-56787-200)

Ear tags: 1 to 100 (EP-56780-100), 1 to 100 (EP-56782-100), 1 to 100 (EP-56783-100), 1 to 100 (EP-56781-100), 1 to 100 (EP-56784-100), 101 to 200 (EP-56780-200), 101 to 200 (EP-56782-200), 101 to 200 (EP-56783-200), 101 to 200 (EP-56781-200), 101 to 200 (EP-56784-200), 1 to 100 (EP-56786-100), 1 to 100 (EP-56788-100), 1 to 100 (EP-56789-100), 1 to 100 (EP-56790-100), 1 to 100 (EP-56787-100), 101 to 200 (EP-56786-200), 101 to 200 (EP-56788-200), 101 to 200 (EP-56789-200), 101 to 200 (EP-56790-200), 101 to 200 (EP-56787-200)

LAB