


## La Pias Ear Tag - 5 colors



La Pias Ear Tag, selection of 5 different colors,  
numbered - box of 100 pce

**Product codes:**

Reference: EP-56779-20

**Product description:**

Specially suited for use with small lab and research animals (mice, rats and other rodents).  
The light weight easy-to-read tags allow for trouble-free identification.

- Ear tags are **numbered from 1-100** (available up to 999)
- **available in 5 different colors : Green, Blue, Red, Light blue, Black**
- also available as Stainless Steel (non-colored)
- Once applied, they cannot be easily removed
- Laser-etched number provides durability
- ear tag are autoclavable (applicator not)

**Packing : one box of 100 pce**

Material

Numbered Tag: Aluminum

Clasp: Stainless Steel

Size

Mouse: 3.5mm O.D. with 3mm Pin

Rat: 3.5mm O.D. with 5mm Pin

**Product features:**

Material: Tag: Aluminium / Claps: Stainless steel

**Available models**

Animal: Mouse (EP-56779-20), Mouse (EP-56779-40), Rat (EP-56785-20), Rat (EP-56785-40)

Ear tags: 1 to 20, box of 100 (EP-56779-20), 21 to 40, box of 100 (EP-56779-40), 1 to 20, box of 100 (EP-56785-20), 21 to 40, box of 100 (EP-56785-40)