

Tailveiner for mice



Easy to use mouse and rat restrainer specifically designed for tail vein injections

Product codes:

Reference: TV

Product description:

The Tailveiner is the first easy-to-use mouse and rat restrainer specifically designed for tail vein injections. There are 3 sizes :

- SM : 2,86 cm ID (10-30g)
- STD : 3,18 cm (15-40g)
- LG : 3,81 cm (20-50g)

Product features:**Available models**

Color: Red (TV-RED-150-STD), Clear (TV-150-STD), Clear (TV-150-SM), Clear (TV-150-LG), Red (TV-RED-150-SM)

Size Tailveiner Mouse: -STD 3,18 cm ID (TV-RED-150-STD), -STD 3,18 cm ID (TV-150-STD),
-SM 2,86 cm ID (TV-150-SM), -LG 3,81 cm ID (TV-150-LG), -SM 2,86 cm ID (TV-RED-150-SM)



LAB