

## TCAT warming blanket



Maintain temperature after surgery

**Product codes:**

Reference: WB-

### Product description:

This blanket uses far **infrared radiation** to penetrate deep into the body.  
DUal use : **Lays flat or wraps** around the animal body

**Reference Animal Dimensions (cm)** WB-1 Mice and small rats 20 x 15 WB-2 Large rats and bigger animals

25 x 25

The warming blankets feature overtemperature protection which prevents the blanket temperature from exceeding 42°C in the event that the primary feedback sensor becomes disconnected.

Please note that **the blanket must be supplied with the controller TCAT-2** to match the overtemperature unit setting.

### Product features:

### Available models

Animal: Mouse (WB-1), Rat (WB-2)

