

Compact stereotaxic frame- single manipulator for 2 mice



Portable system for two mice simultaneously -
wide working area

Product codes:

Reference: BLAB-ST-8517

Product description:

Compact and portable stereotaxic frame for 2 mice, without a U-frame. The frame is thus lighter and facilitates the access of instrument instruments to the body of the animal. The animal is held firmly and accurately with dedicated nose clip and ear bars (delivered with the instrument). The device is supplied with a micro-manipulator with 3 positioning axes with verniers.

Product features:

Dimensions: 400x255 mm

Available models

Frame type : Manual (BLAB-ST-68517), Digital (BLAB-ST-68518)