

# Flexible Disposable Feeding Needle



box of 100 sterile, flexible PTFE Disposable Feeding Needle.

**Product codes:**

Reference: DT-99

**Product description:**

Plastic (PTFE) animal feeding needles (AFN) allow user to precisely dispense liquid to animal while **minimizing the chances of esophageal piercing**.

sterile, non-toxic, and non-pyrogenic, available in a variety of gauges and lengths.

**Product features:**

Units: Each box has 100 sterile needles

**Available models**

Size: 15G-3" 2.8mm ball (DT-9925), 17G-3.3" 2.4mm ball (DT-9926), 18G-1.5" 2mm ball (DT-9928), 20G-1.5" 2mm ball (DT-9920)