


Small Animal Surgical Board with Magnetic base



The New Standard for small animal surgery, procedure, & dissection platforms.

Product codes:

Reference: GSM1

Product description:

Grace GSM1 Surgery Board **with Magnetic Base System**

Non-skid feet for added stability.

Boards are made of 2,5cm thick high density polyethylene (HDPE) with **excellent chemical resistance, easy-to-clean**, solid work surface.

Size: 30cm (L) x 20cm (W)

Includes :

- 1x Magnetic Stainless Steel Warming Plate (17x10cm),
- 1 Space Gel Heating Pad,
- & Suction cup limb restraints

Includes Magnetic Base System that incorporates corner magnets for holding metal accessories and recessed steel bars for use with magnetic accessories.

Does not includes :

- flexible LED light
- Magnifier

both accessories can be supplied as options

Product features:

Available models

