


Small Animal Surgical Board with Magnetic base- LED and Magnifier



The New Standard for small animal surgery, procedure, & dissection platforms.

Product codes:

Reference: GSM2

Product description:

Grace GSM2 Surgery Board **with LED, Magnifier, & Magnetic Base System**

Non-skid feet for added stability.

Boards are made of 2,5cm thick high density polyethylene (HDPE) with **excellent chemical resistance, easy-to-clean**, solid work surface.

Size: 30cm (L) x 24cm (W)

Includes :

- 1x Magnetic Stainless Steel Warming Plate (17x10cm),
- 1 Space Gel Heating Pad,
- & Suction cup limb restraints

- **Two Flexible LED lights** - provides direct illumination with no shadows and no heat. Flex arm, non-obstructive rimless
- **Magnifying glass** (main lens: 2.3X, inset lens: 4.0X). Allows continuous eye contact with object when moving into magnification area.

Includes Magnetic Base System that incorporates corner magnets for holding metal accessories and recessed steel bars for use with magnetic accessories.

Universal input power supply (100 – 240 VAC, 50-60 Hz). US-plug type, supplied with EU-adapter

Product features:

Available models

