


# Esophageal Probe BIO-ESO-1

Esophageal or rectal applications



**Product codes:**

Reference: BIO-ESO-1

## Product description:

The BIO-ESO-1 Esophageal/Rectal Probe for humans and larger animals.  
Flexible **vinyl covered, soft tipped.**

- Sensor Tip Diameter 7/32"
- Accuracy  $\pm 0.1^{\circ}\text{C}$  in the physiological range
- Time Constant: 5.0 secs.
- Maximum Temperature:  $90^{\circ}\text{C}$
- Minimum Temperature  $0^{\circ}\text{C}$
- Lead Length 5 feet
- Sterilization Protocol Gas ETO, Gluteraldehyde, Isoprophyl Alcohol

## Product features:

## Available models

