

Multi-output anesthesia machine



For simultaneous anesthesia (isoflurane or sevoflurane): 1 to 5 rodents

Product codes:

Reference: R550

Product description:

Complete system with vaporizer, regular assembly, inlet/outlet end caps and tube accessory kit.

You can use this versatile system on 1 to 5 small animals among which **rats**, **mice**, hamsters, guinea pigs, rabbits, cats and other within 7kg. Gas flow of induction chamber can be adjusted independently, ranging from 0 to 2.5L/min. Every mask channel can be controlled independently.

You can choose your air supply (oxygen or air) and your anaesthetic product (isoflurane or sevoflurane).

Product features:

Weight: Less than 7.5kg

Dimensions: 326mm (L) × 180mm (W) × 284mm (H)

Temperature distribution: 10-35°C

Animal: 1 to 5 rats, mice, hamsters, guinea pigs, rabbits, cats and other animals within 7kg

Flow rate: 0.1-4L/min (oxygen flowmeter)

Available models

Anaesthetic: Isoflurane (R550IE), Sevoflurane (R550SE)

