

Induction chambers

Transparent for perfect visualization, compatible with our anesthesia stations

**Product codes:**

Reference: V100

Product description:

The induction chamber is the best choice for rapid induction of anesthesia, by inhalation of anesthetic gases, before animal surgery. It takes only 2 to 5 minutes. It is also very useful when you can't use the masks with awake animals. The transparent plexiglass material is convenient to observe in real time the state of animals.

The anesthetic gas is inducted into the chamber, and exhausted to the connected charcoal canister.

Product features:

Material: Plexiglass

Animal: 4 different sizes: mouse, rat, rabbit/cat and large animal

Colors: Transparent

Diameter: 6 mm input ID - 22 mm output ID

Thickness: 1.2 cm (upper cover) - 0.5 cm (other)

Available models

Animal: Rat (V101), Mouse (V100), Rabbit (V102)

