

## Cordless trimmer for rodent BIO-1584



Robust cordless trimmer, offering excellent handling for rodent

**Product codes:**

Reference: BIO-1584

**Product description:**

Cordless trimmer, equipped with a quiet and powerful motor, allowing very fine cuts on a wide variety of animals.

It can be used cordless or with the power cord.

Super light and quiet, only 120g. Ideal for handling rats and mice.

Quick charge : 45 minutes

Delivered with its surgical cutting head **reference TC1584**.

For easy mowing of very small animals and narrow areas, we have an optional precision blade set **reference TC1584P**.

**Product features:**

Width: 32 mm  
Weight: 120 gr  
Charging time: Quick charge : 45 min  
Cutting length: 0.4 mm  
Battery: 75 min

**Available models**

