

# Ear Punch Pincher



Ear punch with slanted edge, for identification and/or sample collection

**Product codes:**

Reference: EP-AT-7000

**Product description:**

This ear punch is used for identification and/or sample collection.

The angled edge allows cutting rather than puncturing, to minimize animal discomfort.

Requires very gentle pressure to cut the hole.

**Reduces musculoskeletal disorders**

Biopsy sample falls perfectly into the attached 1.5 ml tube.

**Autoclavable** but must be lubricated after use for long and trouble-free use.

Made of stainless steel

Sold without tubes, in option we sell the tubes in set of 100 pcs, ref : EP-30120086.

**Product features:**

Diameter: Hole : 2mm or 2,75mm

**Available models**

Ear punch with vial holder: 2 mm punch (EP-AT-7000), 2,75 mm punch (EP-AT-7075)



LAB