

Ear Punch Pincher



Ear punch with slanted edge, for identification and/or sample collection

Product codes:

Reference: EP-AT-7000

Product description:

This ear punch is used for identification and/or sample collection.

The angled edge allows cutting rather than puncturing, to minimize animal discomfort.

Requires very gentle pressure to cut the hole.

Reduces musculoskeletal disorders

Biopsy sample falls perfectly into the attached 1.5 ml tube.

Autoclavable but must be lubricated after use for long and trouble-free use.

Made of stainless steel

Sold without tubes, in option we sell the tubes in set of 100 pcs, ref : EP-30120086.

Product features:

Diameter: Hole : 2mm or 2,75mm

Available models

Ear punch with vial holder: 2 mm punch (EP-AT-7000), 2,75 mm punch (EP-AT-7075)



LAB