

# Gas evacuation apparatus



Allow to safely and actively remove anesthesia waste gas

**Product codes:**

Reference: R546Pro

**Product description:**

When using up to 5 masks, the anesthetic gases are pulled through an activated carbon canister so that the cleaned air can be returned to the room.

It avoids the utilisation of a venting hood while working safely.

**Product features:**

Noise: ?50dB

Dimensions: 215 mm (W) x 215 mm (D) x 170 mm (H)

Other features: 3 levels of overweight alarm: 990, 1010 and 1500 g

Display: Weight of the canister and suction speed

Norms compliance: CE certification

Flow rate: Can be adjusted depending on the stations (masks)

**Available models**

# O LAB