

ThermoStar Homeothermic Monitoring System



ThermoStar Temperature Controller with 1x 12x20,5cm heating pad and 1x 1,5mm diameter temperature probe

Heating range: 20-45°C - 0,1°C resolution

Product codes:

Reference: 6902-

Product description:

Product features:

Weight: 1,3kgs

Dimensions: 205*148*80mm

Power supply: 100-240V 50/60Hz

Resolution: 0,1°C

Working temperature: 20° to 45°C

Available models

Size: S : 7cm x 10cm (69027), M : 9cm x 17cm (69026), L : 12cm x 20,5cm (69020)