

Animal marker with a large tip



Animal marker for **quick and easy** identification with a **standard tip**

Product codes:

Reference: BIO-504-

Product description:

The marker with **large tip** allows **quick and easy** identification of rats and mice

For use on the rodent's coat. Color will not remain on the animal's skin

These pigments are non-toxic and free of any hazardous materials

Possibility to mark approximately 300 animals (depending upon animal size and/or applied area)

Available in different colors: green, blue, red, black, yellow, orange

Tip size : 14.7 x 6.62 mm

Product features:

Available models

Color: Blue (BIO-504-401), Green (BIO-504-411), Red (BIO-504-421), Orange (BIO-504-441), Yellow (BIO-504-431), Black (BIO-504-501)

Packaging: unit (BIO-504-401), unit (BIO-504-411), unit (BIO-504-421), unit (BIO-504-441), unit (BIO-504-431), unit (BIO-504-501)

