

Animal marker with a fine tip



Animal marker for **quick and easy** identification with a **fine tip**

Product codes:

Reference: BIO-504

Product description:

The marker with **fine tip** allows **quick and easy** identification of rats and mice

For use on the rodent's coat. Color will not remain on the animal's skin

These pigments are non-toxic and free of any hazardous materials

Possibility to mark approximately 300 animals (depending upon animal size and/or applied area)

Available in different colors: green, blue, red, black, yellow, orange

Product features:

Available models

Color: Red (BIO-504-471), Black (BIO-504-511), Orange (BIO-504-491), Blue (BIO-504-451),
Green (BIO-504-461), Yellow (BIO-504-481)

Packaging: unit (BIO-504-471), unit (BIO-504-511), unit (BIO-504-491), unit (BIO-504-451), unit
(BIO-504-461), unit (BIO-504-481)

