

Accessories for Rapid Tags

Applicator or Removal Tool for tag Rapid Lab



Product codes:

Reference: BIO-RAPID

Product description:

Applicator or Removal Tool for tag Rapid Lab

Applicator : Stainless Steel and Acetal , all components are autoclave safe at 121 degrees Celsius

Remove Tool : is not intended for autoclave but may be disinfected by alcohol

Product features:

Available models

Rapid: removal tool (BIO-RAPIDREMOV), applicator (BIO-RAPIDAPPLI)