

## BioLite Stand



The perfect stage for intubating mice and rats of all sizes

**Product codes:**

Reference: RIS-200

**Product description:**

The stand is constructed from a single piece of polycarbonate plastic with a seamless design for intubating rats of all sizes.

The top of the stand provides a hand rest for stabilization during intubation that is designed to support significant weight.

The compact seamless design and its construction from durable polycarbonate make this stand easy to clean and will not crack with the use of harsh cleaning agents as is common with similar stands constructed from acrylic. This product provides a simple but durable stand that is easy to clean, store, and move to any location.

Elastics included

**Product features:**

**Available models**

Animal: Rat (RIS-200), Mouse (RIS-100)

