

Red induction chamber for mice and rats

Red chamber, compatible with our anesthesia stations



Product codes:

Reference: V105

Product description:

The induction chamber is the best choice for rapid induction of anesthesia, by inhalation of anesthetic gases, before animal surgery. It takes only 2 to 5 minutes. It is also very useful when you can't use the masks with awake animals.

The anesthetic gas is inducted into the chamber, and exhausted to the connected charcoal canister.

Dim : 31cm*15cm*19cm

Product features:

Material: Plexiglass

Diameter: 6 mm input ID - 22 mm output ID

Thickness: 1.2 cm (upper cover) - 0.5 cm (other)

Available models

O LAB