


Rabbit auricular vein training simulator



Ideal rabbit auricular simulator for learning the skills required for humane and ethical handling and procedural competencies, without the need to use a live animal

Product codes:

Reference: VTQ-RABBIT-EAR

Product description:

BIOSEB follows and offers for many years alternatives to the use of live animals for training in surgical procedures or in an academic and applied work.

Our rabbit auricular training model is a purposely designed training aid to provide the skills required for humane and ethical handling and procedural competencies, without the need to use a live animal.

The soft and flexible rabbit ear has an integrated vein tube, which allows training of injection techniques on the auricular vein in a lifelike manner.

It allows to avoid stress to both animal and handler.

This model is ideal for both veterinary and biomedical establishment

specimen.

Product features:

Available models

